

# BONFIG


## laser cube

Container Application:  
Ampoules, Vials (up to 250 ml),  
Cartridges, Pre-filled syringes  
Products: freeze dried, liquid, powder.



# BONFIG laser cube

Bonfig Laser cube is a benchtop system for performing Headspace Gas Analysis (HGA) of sterile pharmaceutical containers.



It is a non-destructive, non-invasive laser-based inspection technology for measuring the headspace levels of gases, such as oxygen and carbon dioxide, as well as monitoring moisture levels.

The inspection method is created to meet specific customer requirements offering extreme stability and accuracy in inspection even where the headspace is limited. Bonfig Laser Cube is a compact and lightweight system that is easy to use and set up via integrated PC and any wireless touchscreen tablet.

HGA inspection process is based on the Tunable Diode Laser Absorption Spectroscopy (TDLAS) method which uses a laser beam to detect the target molecules within container headspace. HGA is therefore ideal for the accurate investigation of:

1. Translucent containers
2. Headspace conditions for products packaged under modified atmosphere
3. Closure integrity in pharmaceutical finished containers



## Features & Benefits

- Non-invasive and non-intrusive inspection method
- The technology is designed to make Etalon effect eligible
- Fast, reliable and repeatable results
- Nitrogen purging is not required during oxygen measuring
- Compact and durable design
- Maintenance-free & Quick changeover
- Cost-effective solution with low power consumption
- Fully automated test cycle sequencing with manual loading and unloading of containers.
- Data can be stored and exported (production, raw data, events, alarms)
- Enhanced easy-to-use HMI integrated functions.
- HMI real time display of statistics and raw data
- Computerised system is designed to comply with FDA 21 CFR Part 11 and EU Annex 11
- Validation package guarantees complete and efficient regulatory compliance



## Headspace Inspection Capabilities:

Oxygen level  
Moisture level  
Carbon Dioxide level  
Absolute Pressure valve

## Machine Operation

1. Load container into the inspection position
2. Close chamber
3. Push the start button on the HMI to start the test cycle
4. Measurement of target gas within container headspace is taken
5. The HMI displays clearly on the screen whether the container is conforming or non-conforming
  - Containers are classified as non-conforming if the target gas level within the headspace has been exceeded
  - In all other cases containers are classified as conforming.
6. The container is removed and managed in line with internal procedures depending on inspection outcome

## Technical Specification

<b>PACK TYPE:</b>	Ampoules, vials (up to 250 ml), Cartridges, Pre-filled syringes
<b>OPERATION:</b>	Semi-automatic
<b>CONSTRUCTION:</b>	Anodised Aluminium and ABS/Perspex top Cover
<b>USER INTERFACE:</b>	10" MEDICAL TABLET PC
<b>UTILITIES:</b>	<ul style="list-style-type: none"> <li>— Electrical 24v DC 4A</li> <li>— Compressed Air: 6 Bar</li> </ul>
<b>CONFIGURATION:</b>	<ul style="list-style-type: none"> <li>— 3 USB output ports (1 port dedicated to printer or to data download)</li> <li>— 2 Ethernet</li> <li>— 1 HDMI video</li> </ul>
<b>TEST CYCLE:</b>	From 5 seconds
<b>TOOLING CHANGEOVER:</b>	N/A
<b>AUDIT COMPLIANCE:</b>	FDA 21 CFR Part 11 and EU Annex 11
<b>MACHINE DIMENSIONS:</b>	403 mm x 356 mm x 428 mm
<b>WEIGHT:</b>	Machine: 23kg Shipping Weight: 65kg
<b>WARRANTY:</b>	Supplied with a 12 month warranty. (Service Level Agreements and/or extended warranties are available for additional support.



**TASITEST**  
PACKAGING TEST  
& INSPECTION

via Rondona, 33  
44049 - Vigarano Pieve (FE)  
Italy

☎ +39 0532 715 631

2440 W Corporate Preserve Drive,  
Suite 600, Oak Creek,  
WI 53154 USA

☎ +1 4146 713 332

Unit 25 Carrowreagh Business Park  
Carrowreagh Road Dundonald  
Belfast,  
BT16 1QQ  
United Kingdom

☎ +44 2890 484 848

Liebigstrasse 5  
D-85757 Karlsfeld  
Germany

☎ +49 8131 593 910

Dongfang Rd 3601,  
Harvest Industrial Park  
No.7 Building, Room 403,  
Pudong, Shanghai, 200125  
China

☎ +86 2158 366 290