

Protocol for cryopreservation of RPTEC/TERT1 cells

Version: January 2019

Reagents
<p>Trypsin-EDTA (Gibco, Cat# 25300054) Trypsin-Inhibitor (Gibco, Cat# R007100) Freezing medium: CryoStor® cell cryopreservation medium CS10 (Sigma-Aldrich, Cat# C2874) Storage temperature: liquid nitrogen</p> <p>ProxUp ready-to-use medium (Evercyte, Cat# MHT-003) which contains:</p> <ul style="list-style-type: none"> • DMEM/Ham's F-12 (1:1) (Biochrom, Cat# F4815) • 10 mM HEPES-buffer (Biochrom, Cat# L1613) • 2 mM GlutaMAX™-I (Gibco, Cat# 35050-038) • 10 ng/ml hEGF (Sigma-Aldrich, Cat# E9644) • 5 pM 3,3',5-Triiodo-L-thyronine sodium salt (Sigma-Aldrich, Cat# T6397) • 3.5 µg/ml L-Ascorbic Acid (Sigma-Aldrich, Cat# A4544) • 5 µg/ml Transferrin Holo (Merck Millipore, Cat# 616424) • 25 ng/ml Prostaglandine E1 (Sigma-Aldrich, Cat# P8908) • 25 ng/ml Hydrocortisone (Sigma-Aldrich, Cat# H0396) • 8.65 ng/ml Sodium-Selenite (Sigma-Aldrich, Cat# S5261) • 100 µg/ml G418 (InvivoGen, Cat# ant-gn-5) • 5 µg/ml Insulin (Sigma-Aldrich, Cat# I9278) <p>Ethanol, 70 %</p>
Practical application
<p>Freezing of cells: Detach the cells from the culture vessel by using Trypsin and Trypsin-Inhibitor as described above, resuspend the detached cells in growth medium and centrifuge at 170 g for 5 min. Then, discard the supernatant, resuspend in the remaining droplet and add freezing medium (4°C) to reach a cell density of about 1.5 - 2 x 10⁶ cells/ml (for thawing in a 25 cm² culture flask). Transfer 1 ml of this cell suspension to each pre-cooled cryovial and immediately transfer the cells to -80°C. After 24 hours transfer the vials to the liquid nitrogen tank.</p> <p>Thawing of cells: Add 6 ml of growth medium to a 25 cm² culture flask and place the culture flask in the incubator for at least 30 min. Take a vial of frozen cells, rinse outside with ethanol and pre-warm in hand until one last piece of frozen cells is seen. Then, immediately transfer the content of the vial to a 15 ml centrifugation tube pre-filled with 9 ml of medium pre-cooled to 4°C and centrifuge for 5 min at 170 g. Then, discard the supernatant and resuspend the cells in the remaining droplet. Add 1 ml of pre-warmed</p>

medium to the cells, transfer the cell suspension to the prepared culture flask and incubate at 37°C in a suitable incubator.

RPTEC/TERT1 cells are grown in ProxUp medium at 37°C in a humidified atmosphere with 5 % CO₂. Perform a medium change 24 hours after thawing. If the cells are already confluent at this point, they should be passaged (see Evercyte protocol for in vitro propagation of RPTEC/TERT1 cells). For the first passages after thawing we recommend splitting not higher than 1:2.

Related products

- ProxUp ready-to-use medium, 500 ml (Evercyte, Cat# MHT-003)
- ProxUp basal medium, 500 ml (Evercyte, Cat# MHT-003-B)
- ProxUp supplements, 9 vials (Evercyte, Cat# MHT-003-S)