

Protocol for cryopreservation of HCEC-1CT cells

Version: April 2017

Material
<p>Primaria™ culture flasks (Corning, Cat# 353808 and Cat# 353810)</p> <p>Trypsin-EDTA (Gibco, Cat# 25300054) Trypsin-Inhibitor (Gibco, Cat# R007-100) Freezing medium which contains:</p> <ul style="list-style-type: none"> • ColoUp ready-to-use medium (Evercyte, Cat# MHT-039) • 10 % DMSO • 10 % Cosmic Calf Serum (Hyclone, Cat# SH30087) <p>Storage temperature: liquid nitrogen</p> <p>ColoUp ready-to-use medium (Evercyte, Cat# MHT-039) which contains:</p> <ul style="list-style-type: none"> • DMEM / Medium 199 Earle's, 4+1 (Biochrom, Cat# F0435 and Cat# FG0615) • 4 mM GlutaMAX™-1 (100X), (Gibco, Cat# 35050-038) • 2 % Cosmic Calf Serum (Hyclone, Cat# SH30087) • 20 ng/ml EGF (Sigma-Aldrich, Cat# E9644) • 10 µg/ml Insulin (Sigma-Aldrich, Cat# I9278) • 2 µg/ml Apo-Transferrin (Sigma-Aldrich, Cat# T2036) • 5 nM Sodium-Selenite (Sigma-Aldrich, Cat# S5261) • 1 µg/ml Hydrocortisone (Sigma-Aldrich, Cat# H0396) <p>Ethanol, 70%</p>
Practical application
<p>Freezing of cells: Cells are detached from the culture vessel by using Trypsin-EDTA and Trypsin-Inhibitor as described in Evercyte's protocol for in vitro propagation of HCEC-1CT cells and centrifuged at 170 g for 5 min. Then, the supernatant is discarded, the cell pellet is resuspended in the remaining droplet and freezing medium (pre-cooled to 4°C) is added to reach a cell density of about 1 - 2 x 10⁶ cells/ml (for thawing in a 25 cm² Primaria™ culture flask). Then, 1 ml of this cell suspension is added to each pre-cooled cryovial which are then immediately transferred to -80°C. After 24 hours the vials are transferred to liquid nitrogen for long-term storage.</p> <p>Thawing of cells: 6 ml of growth medium are added to a 25 cm² culture flask, which is transferred for at least 30 min to a humidified incubator to allow the medium to reach 37°C and its normal pH. Then, a vial of frozen cells is taken, rinsed outside with ethanol and pre-warmed in the hand until one last piece of frozen cells is seen. Thereafter, the content of the vial is immediately transferred to a 15 ml centrifugation tube pre-filled with 9 ml of medium pre-cooled to 4°C and centrifuged for 5 min at 170 g. Then, the</p>

supernatant is discarded, the cell pellet is resuspended in the remaining droplet and 1 ml of the pre-warmed medium is added to the cells. This cell suspension is then transferred to the prepared culture flask and incubated at 37°C in a suitable incubator. After 24 h a medium change is performed. If the cells are already confluent at this point, they should be passaged (see Evercyte's protocol for in vitro propagation of HCEC-1CT cells).

Related products

- ColoUp ready-to-use medium, 500 ml (Evercyte, Cat# MHT-039)