

Low volume HIV-1 neutralisation assay

Issue

The evaluation of vaccine efficacy in humans during clinical studies calls for assessing the presence of HIV neutralising antibodies in vaccinees. However, *in vitro* testing of precious biological samples available in limited quantity may be challenging. Here we describe an adaptation of our *in vitro* neutralisation assay to accommodate such low sample volumes.

Solution

In our modified format, antibody source (plasma, sera, ...) is first incubated with virus and reporter cells^(*) in a 20 μ L volume for several hours. Afterwards, antibodies and viruses are removed, replaced with normal culture medium and cell culture is resumed for several days before to proceed to reporter quantitation.

^(*): In our reporter cell line ($CD4^+$, $CXCR4^+$, and $CCR5^+$) HIV replication is monitored by a quantitative reporter induction.

Validation

Neutralisation assay was performed using both normal and modified format. The inhibition parameters are reported in the table below:

	High volume		Low volume	
	EC ₅₀	EC ₉₀	EC ₅₀	EC ₉₀
2F5	6	165	10	186
4E10	3.4	233	4.7	228

The tested modification led to non-significant variations, within the assay limits. Therefore, it was further assessed using clonal viruses of different subtypes and plasma from HIV-infected long-term non-progressor patients (see figure 1).




Figure 1. Validation of the modified assay using HIV-1 clones of various subtypes and plasma from long-term non-progressors. Inhibition curves were fitted using Xlfit software and inhibition parameters derived from the fitted curves.

Conclusion

We have successfully adapted our HIV neutralisation assay to accommodate low sample volume. Using this format, 6 μ L of biological samples allow for the testing of 4 recombinant viruses along with 2 specificity controls (dose-response curves).