



SPECIFICATIONS

TAP[®] TARTARIC ACID PELLETS (USP/NF / JP / Ph. EUR)

Type	TAP 200	TAP 300-425	TAP 400	TAP 425-500	TAP 425-600	TAP 500	TAP 600	TAP 700	TAP 850-1000	Test Method
Appearance	White or nearly white or beige particles									Visual
Physical Parameters										
Particle Size	100 – 200 µm	300 – 425 µm	300 – 500 µm	425 – 500 µm	425 – 600 µm	400 – 600 µm	500 – 710 µm	600 – 800 µm	850 – 1000 µm	Ph. EUR
	≥ 75 %									
Bulk Density	0.80 - 1.1 g/ml (for information only)									Ph. EUR
Identification										
A Solution strongly acidic	strongly acidic									Ph. EUR
B Reactions of tartrates	positive									Ph. EUR / USP/NF
C Infrared Absorption	IR-spectrum (ATR) equivalent reference standard									USP/NF
Chemical Parameters										
Appearance of solution	Clear and not more intensely colored than the reference solution Y									Ph. EUR
Water content	≤ 0.5%									Ph. EUR
Loss on drying	≤ 0.5% (for information only)									USP/NF
Calcium	≤ 200 ppm									Ph. EUR
Chloride	≤ 100 ppm									Ph. EUR
Specific optical rotation	+12.0 to 12.8°									Ph. EUR / USP/NF
Limit of oxalate	no turbidity									USP/NF
Residue on ignition	≤ 0.05%									Ph. EUR / USP/NF
Sulfate	no turbidity									USP/NF
Arsenic	≤ 1 ppm									JP
Assay	99.7 to 100.5 %									USP/NF
Microbiological Parameters										
TAMC	< 1000 CFU/g									Ph. EUR / USP/NF / JP
TYMC	< 100 CFU/g									Ph. EUR / USP/NF / JP
E. coli	absence in 1 g sample									Ph. EUR / USP/NF / JP

Organic solvents in accordance to Ph.Eur., 5.4 and USP <467> (CPMP/ICH/283/95) are not used neither by manufacturing of TAP[®] nor by cleaning of equipment. The starting material of TAP[®] is exclusively of vegetable origin. A contamination with animal material by manufacturing, storage or shipment in the original closed containers will not occur. Therefore, all types of TAP[®] are free from Bovine Spongiform Encephalopathy (BSE) and Transmissible Spongiform Encephalopathy (TSE). All types of TAP[®] comply with the requirements of EMA/410/01 rev.3

TAP is a proprietary registered mark belong to third parties and not to Pharmatrans

Contact

Pharmatrans Sanaq AG
Gewerbestrasse 18 - CH-4123 Allschwil
Tel.: +41 61 225 90 00
Fax: +41 61 225 90 01
Email: info@pharmatrans-sanaq.com
Web: www.pharmatrans-sanaq.com

21.01.20