

# Human dermal keratinocytes NHEK/SVTERT3-5

Good experiments start with the right choices – hTERT immortalized cell lines retain the cell-type specific phenotype while constantly growing. No more lot-to-lot variability. No more growth arrest.

Just the perfect choice!



www.evercyte.com

## Human dermal keratinocytes (NHEK/SVTERT3-5)

Human keratinocytes represent the major component of the epidermal tissue with an essential role in forming an effective barrier between the human body and the outside.

## NHEK/SVTERT3-5 in a nutshell

- Original tissue / isolated cells: human adult skin (pendulous abdomen) / keratinocytes (NHEK)
- Life span extension of isolated NHEKs by introduction of hTERT (catalytic subunit of telomerase) and SV40 ER (Simian Virus 40 early region / largeT and small t antigens)
- Single clone with distinguished keratinocyte markers and functions and unlimited growth potential
- · Ability to differentiate into a well-organized 3D skin equivalent in an air-liquid interface
- Expression of typical keratinocyte markers after 3D differentiation


## **Cell-type specific characteristics**

#### Continuous growth in vitro

The cell line was continuously cultured for more than 50 population doublings without showing signs of growth retardation or replicative senescence with a constant population doubling time of 48-60 hours.

#### Marker Expression in 3D skin equivalents

Immunofluorescence stainings of 3D skin equivalents established with NHEK/SVTERT3-5 cells show expression of the cell-type specific markers Keratin-1, Keratin-2, Keratin-10 and Loricrin with a staining pattern similar to normal skin.



## Applications

- Study of pathogenesis of skin-related diseases
- Representative in vitro model to study wound healing processes
- Establishment of standardizable 3D skin equivalents for toxicity studies
- Study of drug delivery across skin barrie skin irritation, corrosion
- Establishment of gene-edited in vitro model systems



## Adherence to GCCP-Standards!

Evercyte is committed to follow the principles of Good Cell Culture Practice (GCCP, Coecke et al., 2005). Therefore, our cell lines are:

- established following highest ethical standards (studies are approved by IRB in accordance with the Declaration of Helsinki)
- **quality tested** (sterility, absence of specific humanpathogenic viruses, STR-Profile, longevity)
- characterized for expression of cell type specific markers and functions

